



Entod Research Cell UK Ltd.
Quality, Innovation & Research

I-DEW FLO

FLUORESCEIN OPHTHALMIC STRIPS FOR CORNEAL STAINING

1mg Fluorescein Strips

The Corneal Staining Expert

- Permeates into the Intercellular Space associated with any Corneal Epithelial Cellular Disruption
- Useful to **Determine the Site of an Ocular Injury**
- **Allows observation** of Superficial Punctate Keratitis, Superior Epithelial Arcuate Lesion, Inferior Arcuate Staining and Three & Nine o'clock staining.
- Useful for **Evaluating Contact Lens Fitting & Comfort and Contact Lens-related Dryness**


Optimum Design & Premium Quality

- 100 Sterile Strips individually wrapped in Transparent pouches
- Produced using the Highest Grade Fluorescein & Premium Paper
- Fully Automated Production with "NO HUMAN TOUCH"
- Designed in the UK & Manufactured under GMP



AS RECOMMENDED BY

tfosDEWS II



Slit Lamp Set-Up:

Wratten #12 filter,
Cobalt Light,
High Illumination

For more Information Contact:

Entod Research Cell UK Ltd.

15 Tottenham Lane, Hornsey, London N8 9DJ, UK.

Tel: (+44) 0208 144 0963

Email: contact@ophthalmicsuk.com