

# I-DEW GREEN



Entod Research Cell UK Ltd.

Quality, Innovation & Research

## LISSAMINE GREEN OPHTHALMIC STRIPS FOR CONJUNCTIVAL STAINING

### 1.5 mg Lissamine Green Strips

CE  
0068  
Certified

AS RECOMMENDED BY


  
tfos DEWS II

### The Conjunctival Staining Expert

- Stains Dead and Degenerate Conjunctival Epithelial Cells, yet does not stain healthy cells
- Same Staining Profile as Rose Bengal but **better tolerated by patients**
- **Allows observation** of Bulbar Staining Circumlimbal Staining, Lid Wiper Epitheliopathy & Lid Parallel Conjunctival Folds
- Useful for **Testing Early Dry Eyes and Contact Lens-related Dryness**
- Can **Evaluate Contact Lens Fitting & Comfort**

### Optimum Design & Premium Quality

- 100 Sterile Strips individually wrapped in Transparent pouches
- Produced using the Highest Grade Lissamine Green & Premium Paper
- Fully Automated Production with "NO HUMAN TOUCH"
- Designed in the UK & Manufactured under GMP



**Slit Lamp Set-Up:**

White Light,  
Low Illumination

**For more Information Contact:**

## Entod Research Cell UK Ltd.

15 Tottenham Lane, Hornsey, London N8 9DJ, UK.

Tel: (+44) 0208 144 0963

Email: [contact@ophthalmicsuk.com](mailto:contact@ophthalmicsuk.com)