

Ophthalmic instruments tables

# MD- series

Slim line



# MD-series

Slim Line



## MD-1

Table-top: 660mm x 430mm  
Max.load: 60 kg



## MD-2

Table-top 860mm x 430mm  
Max.load: 100 kg



## COMBO-1

Table-top: 660mm x 430mm  
Max.load: 60 kg

### Technical specification:

- electrically adjustable table-top height: 610mm - 895mm (+/- 5%)
- input power: 110 - 230V 50 - 60 Hz
- 2 accessory sockets standard
- robust base with 2 casters/2 levelling screws standard (without MD-3 table)
- MD-3 with two paralleled elevating columns
- each column with robust metal base with two inhibitory casters

# MD-series

Slim Line



## COMBO-2

Table-top: 880mm x 430mm  
Max.load: 100 kg



## MD-V

Table-top: V-shape 1040mm x 550mm  
Max.load: 100 kg



## MD-3

Table-top: 1230mm x 550mm  
Maximum load: 150 kg

### Technical specification:

- electrically adjustable table-top height: 610mm - 895mm (+/- 5%)
- input power: 110 - 230V 50 - 60 Hz
- 2 accessory sockets standard
- robust base with 2 casters/2 levelling screws standard (without MD-3 table)
- MD-3 with two paralleled elevating columns
- each column with robust metal base with two inhibitory casters